



ALDOLASE CALIBRATOR

Animal calibration serum for Aldolase determination

INTENDED USE

Concentrations and activities of components are studied for an excellent use even in automatic system of clinical chemistry, even in manual procedures.

PRESERVATION AND STABILITY

Preserve the product not reconstituted at 2-8°C, stable until expiration date on the package.
Calibrator reconstituted is stable 5 days at 2-8°C.

MATERIALS REQUIRED BUT NOT SUPPLIED

Automatic Micropipette, normal laboratory equipment, distilled or deionised water.

REAGENTS PREPARATION

Open carefully one vial of liophile, trying to avoid loss of material, and with the enclosed pipette add exactly 1.0 ml of distilled water, for this is better employing an automatic pipette. Close with attention and let the liophile recovers for 30 minutes, shake carefully trying to avoid foam formation.

PRECAUTIONS

- Product contains Sodium Azide 0.095% as preservative.
- Avoid contact with skin and mucous membrane.
- Use the normal precautions for lab behaviour.
- For in vitro use only.

WASTE DISPOSAL

Product is intended for professional laboratories. Waste products must be handled as per relevant security cards and local regulations.

PACKAGING

CODE CC03340

Calibrator 2 x 1 ml (liophile)

REFERENCES

Archives data.

CALIBRATION VALUES

LOT C108.12

2014/04

COMPONENT	T (°C)	CALIBRATION VALUE	UNIT
CALIBRATOR	37°C	8.99	U/l

MANUFACTURER

LTA s.r.l.
Via Milano 15/F
20060 Bussero (Milan) ITALY
Tel: ++39 02 95409034
Fax: ++39 02 95334185
e-mail: info@ltaonline.it
Website: <http://www.ltaonline.it>

SYMBOLS

	Only for IVD use
	Lot of manufacturing
	Code number
	Storage temperature interval
	Expiration date
	Warning, read enclosed documents
	Read the directions
	Biological risk

Mod. 01.06 (ver. 1.0 - 04/03/2006)

