

CONTROL SET COPPER / ZINC

2 Levels control solutions for determination of Copper and Zinc in biological fluids

INTENDED USE

Concentrations and activities of components are studied for an excellent use even in automatic system of clinical chemistry, even in manual procedures.

PRESERVATION AND STABILITY

Preserve the product at 2-8°C, stable until expiration date on the package.
Once opened controls are stable for 7 days at 2-8°C.

MATERIALS REQUIRED BUT NOT SUPPLIED

Normal laboratory equipment.

REAGENTS PREPARATIONS

Controls are ready to use.
Bring controls at room temperature before use and verify that there isn't crystalline precipitate on the bottom of the bottle. In precipitate presence shake and/or warm a little bit to re-dissolve.

PRECAUTIONS

- The product contains very dilute acids.
- Avoid contact with skin and mucous membrane.
- Use the normal precautions for lab behaviour.
- For in vitro use only.

WASTE DISPOSAL

Product is intended for professional laboratories. Waste products must be handled as per relevant security cards and local regulations.

PACKAGING

CODE CC02780

Control Level 1 2 x 5 ml (liquid)
Control Level 2 2 x 5 ml (liquid)



REFERENCES

Archives data.

CALIBRATION VALUES

LOT K022.18






 2022/01

COMPONENT / METHOD		VALUE	RANGE	UNIT
CONTROL LEVEL 1		Colorimetric tests for the determination of Copper and Zinc in biological fluids		
LOT F018.18  2022/01	COPPER	129	106 – 152	µg/dl
	ZINC	86,8	71,2 – 102,4	µg/dl
CONTROL LEVEL 2		Colorimetric tests for the determination of Copper and Zinc in biological fluids		
LOT G018.18  2022/01	COPPER	335	275 – 395	µg/dl
	ZINC	49,6	40,7 – 58,6	µg/dl

MANUFACTURER

LTA s.r.l.
Via Milano 15/F
20060 Bussero (Milan) - ITALY
Tel: +39 02 95409034
Fax: +39 02 95334185
e-mail: info@LTAonline.it
Website: <http://www.LTAonline.it>

SYMBOLS

- IVD** Only for IVD use
LOT Lot of manufacturing
REF Code number
 Storage temperature interval
 Expiration date (year – month)
 Warning, read enclosed documents
 Read the directions
 Biological risk

Mod. 01.06 (ver. 1.0 – 07/03/2017)

